


other
airway
devices
[created
by Paul
Young
22/12/07]

larynx tube



- Sizes
- Size 3, for Child and Small Adult (Yellow Connector)
 - Size 4, for Medium Adult (Red Connector)
 - Size 5, for Large Adult (Violet Connector)

insertion:



combitube



- features:
- white normal ETT
 - Blue opening distal of the large cuff .
 - Indicators for intubation depth.
 - Cuff blue: 100ml air, seals the hypopharynx.
 - Cuff white: 15 ml air, seals trachea / or esophagus
 - can be inserted blindly
- disadvantages:
- (i) Mostly adult sizes
 - (ii) Potential for esophageal trauma
 - (iii) Problems maintaining seal in some patients

general:

- Cuffed oropharyngeal airway
- Easy to use (like guedel)
- no protection against aspiration
- No protection against gastric inflation

COPA
tube

